



SPRE rabbit pAb

Catalog No	YP-Ab-09501
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPR
Immunogen	Synthesized peptide derived from human SPRE AA range: 148-198
Specificity	This antibody detects endogenous levels of human SPRE
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C79158;FLJ21897;FLJ31917;MGC163164;Spred 2;Sprouty protein with EVH 1 domain 2 related sequence;sprouty related EVH1 domain containing 2;Sprouty related protein with EVH 1 domain 2;zgc:77284;SPRE;SPR;Sepiapterin reductase (SPR) (EC 1.1.1.153)
Observed Band	
Calculated Molecular Weight	47kDa; 48kD
Cell Pathway	Cytoplasm.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
