



REEP6 rabbit pAb

Catalog No	YP-Ab-09502
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	REEP6 C19orf32 DP1L1
Immunogen	Synthesized peptide derived from human REEP6 AA range: 110-160
Specificity	This antibody detects endogenous levels of human REEP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	6111M24Rik;C19orf32;Deleted in polyposis 1 like 1;Dp1l1;FLJ25383; OTTMUSP265;OTTMUSP2651;OTTMUSP2652;OTTMUSP21713; OTTMUSP21714;Polyposis locus protein 1 like 1;Polyposis locus protein 1-like 1; Receptor accessory protein 6;Receptor expression enhancing protein 6; Receptor expression-enhancing protein 6;REEP6;REEP6_HUMAN;TB2 protein like 1;TB2L1;REEP6 C19orf32 DP1L1;Receptor expression-enhancing protein 6 (Polyposis locus protein 1-like 1)
Observed Band	
Calculated Molecular Weight	22 kDa; 23kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, clathrin-coated vesicle membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in circumvallate papillae and testis (PubMed:16720576). Expressed in the retina. Isoform 1 is predominantly present in mature optic cups. Isoform 1 expression is confined to the cell body and inner segment of developing rod photoreceptor cells (PubMed:27889058).
Function	
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images