



RHBL2 rabbit pAb

Catalog No	YP-Ab-09504
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RHBDL2
Immunogen	Synthesized peptide derived from human RHBL2 AA range: 173-223
Specificity	This antibody detects endogenous levels of human RHBL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RHBL2;RHBDL2;Rhomboid-related protein 2 (RRP2) (EC 3.4.21.105) (Rhomboid-like protein 2) [Cleaved into: Rhomboid-related protein 2;N-terminal fragment (NTF);Rhomboid-related protein 2;C-terminal fragment (CTF)]
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Rhomboid-related protein 2, C-terminal fragment]: Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

