



ZGPAT rabbit pAb

Catalog No	YP-Ab-09509
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZGPAT GPATC6 GPATCH6 KIAA1847 ZC3H9 ZC3HDC9 ZIP
Immunogen	Synthesized peptide derived from human ZGPAT AA range: 34-84
Specificity	This antibody detects endogenous levels of human ZGPAT
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	G patch domain-containing protein 6;GPATC6;GPATCH6;KIAA1847;ZC3H9;ZGPAT;ZGPAT GPATC6 GPATCH6 KIAA1847 ZC3H9 ZC3HDC9 ZIP;Zinc finger CCCH-type with G patch domain-containing protein (G patch domain-containing protein 6) (Zinc finger CCCH domain-containing protein 9) (Zinc finger and G patch domain-containing protein)
Observed Band	65 kDa
Calculated Molecular Weight	531 aa, 57 kDa
Cell Pathway	Nucleus .; [Isoform 4]: Nucleus .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
