



ZNF57 rabbit pAb

Catalog No	YP-Ab-09510
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF57 ZNF424
Immunogen	Synthesized peptide derived from human ZNF57 AA range: 246-296
Specificity	This antibody detects endogenous levels of human ZNF57
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger protein 575;ZN575_HUMAN;Znf575;ZNF57;ZNF57 ZNF424;Zinc finger protein 57 (Zinc finger protein 424)
Observed Band	
Calculated Molecular Weight	26 KD; 27kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
