



SUSD4 rabbit pAb

Catalog No	YP-Ab-09514
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SUSD4 UNQ196/PRO222
Immunogen	Synthesized peptide derived from human SUSD4 AA range: 122-172
Specificity	This antibody detects endogenous levels of human SUSD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRO222;sushi domain containing 4;Sushi domain-containing protein 4;UNQ196/PRO222;SUSD4;SUSD4 UNQ196/PRO222
Observed Band	54 kDa
Calculated Molecular Weight	490 aa, 54 kDa
Cell Pathway	[Isoform 1]: Membrane ; Single-pass type I membrane protein .; [Isoform 3]: Secreted .
Tissue Specificity	Isoform 3 is the predominant isoform in all tissues except cortex, cerebellum, kidney, and breast. Isoform 1 is found primarily in the esophagus and the brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

