



THUM2 rabbit pAb

Catalog No	YP-Ab-09515
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	THUMPD2 C2orf8
Immunogen	Synthesized peptide derived from human THUM2 AA range: 299-349
Specificity	This antibody detects endogenous levels of human THUM2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	THUM2;THUMPD2 C2orf8;THUMP domain-containing protein 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in a variety of tissues including brain, colon, gingiva, heart, kidney, liver, lung, placenta, small intestine, spleen and thymus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
