



S29P2 rabbit pAb

Catalog No	YP-Ab-09519
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SNX29P2 RUNDC2C
Immunogen	Synthesized peptide derived from human S29P2 AA range: 105-155
Specificity	This antibody detects endogenous levels of human S29P2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S29P2;SNX29P2 RUNDC2C;Putative protein SNX29P2 (RUN domain-containing protein 2C) (Sorting nexin 29 protein pseudogene 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
