



ERF3B rabbit pAb

Catalog No	YP-Ab-09520
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GSPT2 ERF3B
Immunogen	Synthesized peptide derived from human ERF3B AA range: 182-232
Specificity	This antibody detects endogenous levels of human ERF3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GST2;G1 to S phase transition 2;Eukaryotic peptide chain release factor GTP-binding subunit ERF3B;eRF3b;EC:3.6.5.-;GSPT2 ERF3B;Eukaryotic peptide chain release factor GTP-binding subunit ERF3B (Eukaryotic peptide chain release factor subunit 3b) (eRF3b) (G1 to S phase transition protein 2 homolog)
Observed Band	69-75 kDa
Calculated Molecular Weight	628 aa, 69 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Highly expressed in IUCC stage II colorectal cancer (CRC).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
