



TTC9C rabbit pAb

Catalog No	YP-Ab-09522
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTC9C
Immunogen	Synthesized peptide derived from human TTC9C AA range: 61-111
Specificity	This antibody detects endogenous levels of human TTC9C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Tetratricopeptide repeat protein 9C;TPR repeat protein 9C;TTC9C; Tetratricopeptide repeat protein 9C (TPR repeat protein 9C)
Observed Band	
Calculated Molecular Weight	22 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
