



HEMGN rabbit pAb

Catalog No	YP-Ab-09527
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HEMGN EDAG NDR PRO1037 PRO1620
Immunogen	Synthesized peptide derived from human HEMGN AA range: 419-469
Specificity	This antibody detects endogenous levels of human HEMGN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EDAG;EDAG 1;HEMGN;hemogen;Hemopoietic gene protein;NDR;PRO1620; HEMGN EDAG NDR PRO1037 PRO1620;Hemogen (Erythroid differentiation-associated gene protein) (EDAG-1) (Hemopoietic gene protein) (Negative differentiation regulator protein)
Observed Band	60 kDa
Calculated Molecular Weight	53kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in hematopoietic precursor cells, thyroid and spermatids (at protein level). Expressed in bone marrow, testis, thymus. Expressed in prostate cancer and ovarian cancer. Also expressed in thymus and thyroid tumors, non-Hodgkin lymphoma, various leukemia cell lines, peripheral blood mononuclear cells (PBMCs) and bone marrow mononuclear cells (BMMCs) of patients with leukemia.
Function	
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images