



CTSR2 rabbit pAb

Catalog No	YP-Ab-09529
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CATSPER2
Immunogen	Synthesized peptide derived from human CTSR2 AA range: 189-239
Specificity	This antibody detects endogenous levels of human CTSR2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTSR2;CATSPER2;Cation channel sperm-associated protein 2 (CatSper2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, cilium, flagellum membrane ; Multi-pass membrane protein .
Tissue Specificity	Testis-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
