



SEZ6 rabbit pAb

Catalog No	YP-Ab-09530
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SEZ6
Immunogen	Synthesized peptide derived from human SEZ6 AA range: 108-158
Specificity	This antibody detects endogenous levels of human SEZ6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Brain-specific receptor-like protein C;BSRP-C;hSEZ-6;Seizure protein 6;Seizure protein 6 homolog;Seizure related gene 6;SEZ-6;Sez6;SEZ6_HUMAN;Seizure protein 6 homolog (SEZ-6) (hSEZ-6)
Observed Band	
Calculated Molecular Weight	17kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Localized on dendrites and in the synaptic and postsynaptic fraction. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

