



DK5 rabbit pAbN

Catalog No	YP-Ab-09531
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NME5
Immunogen	Synthesized peptide derived from human NDK5 AA range: 25-75
Specificity	This antibody detects endogenous levels of human NDK5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Inhibitor of p53 induced apoptosis beta;Inhibitor of p53-induced apoptosis-beta; IPIA beta;IPIA-beta;NDK-H 5;NDK5_HUMAN;NDP kinase homolog 5;NM23 H5; nm23-H5;Nm23-M5;NM23H5;NME/NM23 family member 5;NME5;Non-metastatic cells 5;protein expressed in (nucleoside-diphosphate kinase); Nonmetastatic cells 5;protein expressed in;Nonmetastatic protein 23;homolog 5; Nucleoside diphosphate kinase homolog 5;Radial spoke 23 homolog;RSPH23; Testis specific nm23 homolog;Testis-specific nm23 homolog;NDK5;Nucleoside diphosphate kinase homolog 5 (NDK-H 5) (NDP kinase homolog 5) (Inhibitor of p53-induced apoptosis-beta) (IPIA-beta) (Testis-specific nm23 homolog) (nm23-H5)
Observed Band	
Calculated Molecular Weight	24 kDa; 24kD
Cell Pathway	Cell projection, cilium .
Tissue Specificity	Specifically expressed in testis germinal cells.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
