



TSN5 rabbit pAb

Catalog No	YP-Ab-09532
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TSPAN5 TM4SF9
Immunogen	Synthesized peptide derived from human TSN5 AA range: 113-163
Specificity	This antibody detects endogenous levels of human TSN5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TSN5;TSPAN5 TM4SF9;Tetraspanin-5 (Tspan-5) (Tetraspan NET-4) (Transmembrane 4 superfamily member 9)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

