



PPL13 rabbit pAb

Catalog No	YP-Ab-09533
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LGALS14 PPL13
Immunogen	Synthesized peptide derived from human PPL13 AA range: 18-68
Specificity	This antibody detects endogenous levels of human PPL13
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PPL13;LGALS14 PPL13;Placental protein 13-like (Charcot-Leyden crystal protein 2) (CLC2) (Galectin-14) (Gal-14)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
