



AT8B2 rabbit pAb

Catalog No	YP-Ab-09534
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ATP8B2 ATPID KIAA1137
Immunogen	Synthesized peptide derived from human AT8B2 AA range: 908-958
Specificity	This antibody detects endogenous levels of human AT8B2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AT8B2;ATP8B2 ATPID KIAA1137;Probable phospholipid-transporting ATPase ID (EC 3.6.3.1) (ATPase class I type 8B member 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein . Efficient exit from the endoplasmic reticulum requires the presence of TMEM30A or TMEM30B. .
Tissue Specificity	Isoform 3 is ubiquitous, with highest expression in aorta, cerebellum and uterus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
