



RBG1L rabbit pAb

Catalog No	YP-Ab-09536
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RABGAP1L HHL KIAA0471
Immunogen	Synthesized peptide derived from human RBG1L AA range: 639-689
Specificity	This antibody detects endogenous levels of human RBG1L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RBG1L;RABGAP1L HHL KIAA0471;Rab GTPase-activating protein 1-like
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle . Early endosome . Golgi apparatus . Colocalizes on endosomes partially with EEA1 (PubMed:16923123). Colocalizes and cotransports on motile vesicles with ANK2 (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
