



CSMD3 rabbit pAb

Catalog No	YP-Ab-09537
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CSMD3 KIAA1894
Immunogen	Synthesized peptide derived from human CSMD3 AA range: 1623-1673
Specificity	This antibody detects endogenous levels of human CSMD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CSMD3;CSMD3 KIAA1894;CUB and sushi domain-containing protein 3 (CUB and sushi multiple domains protein 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Weakly expressed in most tissues, except in brain. Expressed at intermediate level in brain, including cerebellum, substantia nigra, thalamus, spinal cord, hippocampus and fetal brain. Also expressed in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
