



UN13C rabbit pAb

Catalog No	YP-Ab-09539
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UNC13C
Immunogen	Synthesized peptide derived from human UN13C AA range: 1996-2046
Specificity	This antibody detects endogenous levels of human UN13C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UN13C;UNC13C;Protein unc-13 homolog C (Munc13-3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein . Cell junction, synapse, presynaptic cell membrane ; Peripheral membrane protein . Localized to presynaptic structures. .
Tissue Specificity	Exclusively expressed in brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
