



PEO1 rabbit pAb

Catalog No	YP-Ab-09540
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PEO1 C10orf2
Immunogen	Synthesized peptide derived from human PEO1 AA range: 173-223
Specificity	This antibody detects endogenous levels of human PEO1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATXN8;C10orf2;IOSCA;PEO;PEO1;PEOA3;SANDO;SCA8;Twinkle;Twinkle protein;mitochondrial;TWINL;PEO1 C10orf2;mitochondrial (EC 3.6.4.12) (Progressive external ophthalmoplegia 1 protein) (T7 gp4-like protein with intramitochondrial nucleoid localization) (T7-like mitochondrial DNA helicase)
Observed Band	60 kDa
Calculated Molecular Weight	684 aa, 77 kDa
Cell Pathway	Mitochondrion matrix, mitochondrion nucleoid . Colocalizes with mtDNA in mitochondrial nucleoids, a nucleoproteins complex consisting of a number of copies of proteins associated with mtDNA, probably involved in mtDNA maintenance and expression. .
Tissue Specificity	High relative levels in skeletal muscle, testis and pancreas. Lower levels of expression in the heart, brain, placenta, lung, liver, kidney, spleen, thymus, prostate, ovary, small intestine, colon and leukocytes. Expression is coregulated with MRPL43.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images