



VKIND rabbit pAb

Catalog No	YP-Ab-09541
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KNDC1 C10orf23 KIAA1768 RASGEF2 VKIND hucep-9
Immunogen	Synthesized peptide derived from human VKIND AA range: 1464-1514
Specificity	This antibody detects endogenous levels of human VKIND
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VKIND;KNDC1 C10orf23 KIAA1768 RASGEF2 VKIND hucep-9;Protein very KIND (Cerebral protein 9) (KIND domain-containing protein 1) (Kinase non-catalytic C-lobe domain-containing protein 1) (Ras-GEF domain-containing family member 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, dendrite . Perikaryon .
Tissue Specificity	Expressed specifically in the cerebral cortex.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
