



AbRDH13 rabbit p

Catalog No	YP-Ab-09542
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RDH13 PSEC0082 UNQ736/PRO1430
Immunogen	Synthesized peptide derived from human RDH13 AA range: 1-51
Specificity	This antibody detects endogenous levels of human RDH13
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.1.1.300;PSEC0082;Retinol dehydrogenase 13;SDR7C3;Short chain dehydrogenase/reductase family 7C member 3;RDH13;RDH13 PSEC0082 UNQ736/PRO1430;Retinol dehydrogenase 13 (EC 1.1.1.-)
Observed Band	36 kDa
Calculated Molecular Weight	331 aa, 36 kDa
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein . Localized on the outer side of the inner mitochondrial membrane. .
Tissue Specificity	Widely expressed (PubMed:18039331). In the retina, detected in the inner segment of the photoreceptor cells. Weak signals are observed in a small population of inner nuclear neurons and the inner plexiform layer (PubMed:12226107).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images