



# RIAD1 rabbit pAb

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-Ab-09545                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Reactivity                  | Human;Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Applications                | WB                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Gene Name                   | RIAD1 C1orf230                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Immunogen                   | Synthesized peptide derived from human RIAD1 AA range: 2-52                                                                                                                                                                                                                                                                                                                                                                                   |
| Specificity                 | This antibody detects endogenous levels of human RIAD1                                                                                                                                                                                                                                                                                                                                                                                        |
| Formulation                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                                                       |
| Source                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Purification                | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                                                                                                                                                                                                                                                     |
| Dilution                    | WB 1:1000-2000                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Concentration               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Purity                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Storage Stability           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Synonyms                    | HT5;all trans retinoic acid inducible RING finger;Ari 2;Ari-2;ARI2;ARI2_HUMAN; Ariadne 2;Ariadne 2 protein homolog;ariadne homolog 2 (Drosophila);ariadne homolog 2;Ariadne RBR E3 ubiquitin protein ligase 2;ARIH 2;Arih2;E3 ubiquitin protein ligase ARIH2;E3 ubiquitin-protein ligase ARIH2;protein ariadne 2 homolog;Protein ariadne-2 homolog;Triad 1;TRIAD1;Triad1 protein;UIP48; RIAD1;RIAD1 C1orf230;RIIa domain-containing protein 1 |
| Observed Band               |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Calculated Molecular Weight | 59 kDa; 58kD                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Cell Pathway                |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Tissue Specificity          |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Function                    |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Background                  |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| matters needing attention   | Avoid repeated freezing and thawing!                                                                                                                                                                                                                                                                                                                                                                                                          |

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

---



## Products Images

|  |
|--|
|  |
|  |
|  |
|  |
|  |
|  |
|  |