



BTBD8 rabbit pAb

Catalog No	YP-Ab-09546
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	BTBD8
Immunogen	Synthesized peptide derived from human BTBD8 AA range: 224-274
Specificity	This antibody detects endogenous levels of human BTBD8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BTB (POZ) domain containing 8;BTBD8;BTB/POZ domain-containing protein 8
Observed Band	43 kDa, 34 kDa
Calculated Molecular Weight	278 aa, 43 kDa
Cell Pathway	Nucleus . Localized to nucleus in fetal cells.
Tissue Specificity	Highly expressed in fetal brain. Weakly expressed in adult brain and prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
