



MUC12 rabbit pAb

Catalog No	YP-Ab-09549
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MUC12 MUC11
Immunogen	Synthesized peptide derived from human MUC12 AA range: 847-897
Specificity	This antibody detects endogenous levels of human MUC12
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MUC11;MUC-11;MUC-12;Mucin-11;Mucin-12;MUC12;MUC12 MUC11;Mucin-12 (MUC-12) (Mucin-11) (MUC-11)
Observed Band	550 kDa
Calculated Molecular Weight	558 kDa
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Ubiquitous, with higher expression in colon. Down-regulated in colorectal cancer as well as in the colon of patients with ulcerative colitis (UC) and Crohn's disease (CD).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
