



TRNK1 rabbit pAb

Catalog No	YP-Ab-09550
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TRANK1 KIAA0342 LBA1
Immunogen	Synthesized peptide derived from human TRNK1 AA range: 2101-2151
Specificity	This antibody detects endogenous levels of human TRNK1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRNK1;TRANK1 KIAA0342 LBA1;TPR and ankyrin repeat-containing protein 1 (Lupus brain antigen 1 homolog)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

