



ACO13 rabbit pAb

Catalog No	YP-Ab-09552
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ACOT13 THEM2 HT012 PNAS-27
Immunogen	Synthesized peptide derived from human ACO13 AA range: 47-97
Specificity	This antibody detects endogenous levels of human ACO13
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACO13;ACOT13 THEM2 HT012 PNAS-27;Acyl-coenzyme A thioesterase 13 (Acyl-CoA thioesterase 13) (EC 3.1.2.-) (Thioesterase superfamily member 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol . Mitochondrion . Nucleus . Cytoplasm, cytoskeleton, spindle . During interphase, found both in the nucleus and in the cytoplasm. At mitosis, localizes to the spindle. Colocalizes with tubulin. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
