



TRI55 rabbit pAb

Catalog No	YP-Ab-09553
Isotype	IgG
Reactivity	Human;Rat
Applications	WB
Gene Name	TRIM55 MURF2 RNF29
Immunogen	Synthesized peptide derived from human TRI55 AA range: 75-125
Specificity	This antibody detects endogenous levels of human TRI55
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRI55;TRIM55 MURF2 RNF29;Tripartite motif-containing protein 55 (Muscle-specific RING finger protein 2) (MuRF-2) (MuRF2) (RING finger protein 29)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Nuclear under atrophic conditions and upon mechanical signals. Localizes to the sarcomeric M-band in cardiomyocytes. Colocalizes in part with microtubules (By similarity). .
Tissue Specificity	Highly expressed in muscle. Low-level expression in liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
