



CQ10A rabbit pAb

Catalog No	YP-Ab-09554
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	COQ10A UNQ6192/PRO20219
Immunogen	Synthesized peptide derived from human CQ10A AA range: 136-186
Specificity	This antibody detects endogenous levels of human CQ10A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CQ10A;COQ10A UNQ6192/PRO20219;Coenzyme Q-binding protein COQ10 homolog A;mitochondrial
Calculated Molecular Weight	27kD
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein ; Matrix side .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

