



WHAMM rabbit pAb

Catalog No	YP-Ab-09555
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	WHAMM KIAA1971 WHDC1
Immunogen	Synthesized peptide derived from human WHAMM AA range: 387-437
Specificity	This antibody detects endogenous levels of human WHAMM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1971;WHDC1;WHAMM;WHAMM KIAA1971 WHDC1;WASP homolog-associated protein with actin;membranes and microtubules (WAS protein homology region 2 domain-containing protein 1) (WH2 domain-containing protein 1)
Observed Band	90-100 kDa
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Endoplasmic reticulum-Golgi intermediate compartment . Cytoplasmic vesicle membrane . Golgi apparatus, cis-Golgi network . Localized to a perinuclear compartment near the microtubule-organizing center (MTOC). Also detected on tubulo-vesicular structures in the cell periphery that frequently localized along microtubules. .
Tissue Specificity	Expressed in brain, lung, heart, colon and kidney (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images