



SNX8 rabbit pAb

Catalog No	YP-Ab-09556
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SNX8
Immunogen	Synthesized peptide derived from human SNX8 AA range: 373-423
Specificity	This antibody detects endogenous levels of human SNX8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Mvp1;SNX8;sorting nexin 8;Sorting nexin-8
Observed Band	53 kDa
Calculated Molecular Weight	51kD
Cell Pathway	Early endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Colocalizes with retromer components.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
