



ZN823 rabbit pAb

Catalog No	YP-Ab-09557
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF823 ZFP36
Immunogen	Synthesized peptide derived from human ZN823 AA range: 372-422
Specificity	This antibody detects endogenous levels of human ZN823
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN823;ZNF823 ZFP36;Zinc finger protein 823 (Zinc finger protein ZFP-36)
Observed Band	67kD
Calculated Molecular Weight	67kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
