



HSFY1 rabbit pAb

Catalog No	YP-Ab-09561
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HSFY1 HSF2L HSFY; HSFY2 HSF2L HSFY
Immunogen	Synthesized peptide derived from human HSFY1 AA range: 177-227
Specificity	This antibody detects endogenous levels of human HSFY1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSF2 like;HSF2L;HSFY;HSFY1;HSFY2;HSFY1 HSF2L HSFY;HSFY2 HSF2L HSFY;Heat shock transcription factor;Y-linked (Heat shock transcription factor 2-like protein) (HSF2-like)
Observed Band	45 kDa
Calculated Molecular Weight	401 aa, 45 kDa
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Testis-specific. Present in Sertoli cells and spermatogenic cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
