



ADA12 rabbit pAb

Catalog No	YP-Ab-09563
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ADAM12 MLTN UNQ346/PRO545
Immunogen	Synthesized peptide derived from human ADA12 AA range: 100-180
Specificity	This antibody detects endogenous levels of human ADA12
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A disintegrin and metalloproteinase domain 12;ADA12_HUMAN;ADAM 12; ADAM metalloproteinase domain 12;ADAM12;CAR1;Disintegrin and metalloproteinase domain-containing protein 12;MCMP;MCMPMLtna;Meltrin alpha;Meltrin-alpha;Metalloprotease disintegrin 12 transmembrane;MLTN; MLTNA;OTTHUMP46766;ADA12;ADAM12 MLTN UNQ346/PRO545;Disintegrin and metalloproteinase domain-containing protein 12 (ADAM 12) (EC 3.4.24.-) (Meltrin-alpha)
Observed Band	
Calculated Molecular Weight	97 kDa,68 kDa; 1kD
Cell Pathway	[Isoform 1]: Cell membrane; Single-pass type I membrane protein.; [Isoform 2]: Secreted.; [Isoform 3]: Secreted .; [Isoform 4]: Secreted .
Tissue Specificity	Isoform 1 is expressed in placenta and skeletal, cardiac, and smooth muscle. Isoform 2 seems to be expressed only in placenta or in embryo and fetus. Both forms were expressed in some tumor cells lines. Not detected in brain, lung, liver, kidney or pancreas.
Function	
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images