



NUMBL rabbit pAb

Catalog No	YP-Ab-09564
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NUMBL
Immunogen	Synthesized peptide derived from human NUMBL AA range: 290-340
Specificity	This antibody detects endogenous levels of human NUMBL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Numb-related protein;Numb-like protein;NUMBLIKE;NUMBL;Numb-like protein (Numb-related protein) (Numb-R)
Observed Band	65-70 kDa
Calculated Molecular Weight	65 kDa
Cell Pathway	Cytoplasm . Symmetrically distributed throughout the cytoplasm in non dividing neuroblasts of the CNS. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

