



JAGN1 rabbit pAb

Catalog No	YP-Ab-09566
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	JAGN1
Immunogen	Synthesized peptide derived from human JAGN1 AA range: 111-161
Specificity	This antibody detects endogenous levels of human JAGN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	JAGN1;Protein jagunal homolog 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
