



NTP5 rabbit pAb

C

Catalog No	YP-Ab-09567
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CNTNAP5 CASPR5
Immunogen	Synthesized peptide derived from human CNTP5 AA range: 765-815
Specificity	This antibody detects endogenous levels of human CNTP5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CNTP5;CNTNAP5 CASPR5;Contactin-associated protein-like 5 (Cell recognition molecule Caspr5)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
