



RGPS2 rabbit pAb

Catalog No	YP-Ab-09571
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RALGPS2
Immunogen	Synthesized peptide derived from human RGPS2 AA range: 336-386
Specificity	This antibody detects endogenous levels of human RGPS2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RGPS2;RALGPS2;Ras-specific guanine nucleotide-releasing factor RalGPS2 (Ral GEF with PH domain and SH3-binding motif 2) (RalA exchange factor RalGPS2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cell membrane . Associates with membranes through the PH domain. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
