



APC16 rabbit pAb

Catalog No	YP-Ab-09573
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ANAPC16 C10orf104
Immunogen	Synthesized peptide derived from human APC16 AA range: 16-66
Specificity	This antibody detects endogenous levels of human APC16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANAPC16;Anaphase-promoting complex subunit 16;APC16;CENP-27;Cyclosome subunit 16;ANAPC16 C10orf104;Anaphase-promoting complex subunit 16 (APC16) (Cyclosome subunit 16)
Observed Band	12 kDa
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Chromosome, centromere, kinetochore .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
