



UPP rabbit pAb

Catalog No	YP-Ab-09574
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	UPRT
Immunogen	Synthesized peptide derived from human UPP AA range: 102-152
Specificity	This antibody detects endogenous levels of human UPP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:2.4.2.3;UDRPASE;UP;UPASE;UPase 1;UPP;UPRT;Uracil phosphoribosyltransferase homolog
Observed Band	34 kDa
Calculated Molecular Weight	34 kDa
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Highly expressed in leukocytes, liver, spleen and thymus, with lower expression in brain, lung and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
