



ASXL1 rabbit pAb

Catalog No	YP-Ab-09577
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ASXL1 KIAA0978
Immunogen	Synthesized peptide derived from human ASXL1 AA range: 431-481
Specificity	This antibody detects endogenous levels of human ASXL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ASXL1;KIAA0978;ASXL1 KIAA0978;Putative Polycomb group protein ASXL1 (Additional sex combs-like protein 1)
Observed Band	~200 kDa
Calculated Molecular Weight	84aa,10 kDa; 1541aa,165
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed at low level. Expressed in heart, brain, skeletal muscle, placenta, pancreas, spleen, prostate, small intestine, colon, peripheral blood, leukocytes, bone marrow and fetal liver. Highly expressed in testes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
