



LRC45 rabbit pAb

Catalog No	YP-Ab-09578
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LRRC45
Immunogen	Synthesized peptide derived from human LRC45 AA range: 483-533
Specificity	This antibody detects endogenous levels of human LRC45
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRC45;Leucine-rich repeat-containing protein 45;LRRC45
Observed Band	70 kDa
Calculated Molecular Weight	76 kDa
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Localizes to the proteinaceous linker between the proximal ends of the centrioles. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

