



RN152 rabbit pAb

Catalog No	YP-Ab-09579
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RNF152
Immunogen	Synthesized peptide derived from human RN152 AA range: 150-200
Specificity	This antibody detects endogenous levels of human RN152
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RN152;RNF152;E3 ubiquitin-protein ligase RNF152 (EC 6.3.2.-) (RING finger protein 152)
Observed Band	
Calculated Molecular Weight	22kD
Cell Pathway	Lysosome membrane ; Single-pass membrane protein .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
