



VAMP5 rabbit pAb

Catalog No	YP-Ab-09580
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VAMP5 HSPC191
Immunogen	Synthesized peptide derived from human VAMP5 AA range: 2-52
Specificity	This antibody detects endogenous levels of human VAMP5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Myobrevin;VAMP 5;VAMP5;VAMP5 HSPC191;Vesicle-associated membrane protein 5 (VAMP-5) (Myobrevin)
Observed Band	13-15 kDa
Calculated Molecular Weight	12.8 kDa
Cell Pathway	Cell membrane ; Single-pass type IV membrane protein . Endomembrane system ; Single-pass type IV membrane protein . Golgi apparatus, trans-Golgi network membrane ; Single-pass type IV membrane protein . Associated with the plasma membrane as well as intracellular perinuclear and peripheral vesicular structures of myotubes. Associated with the trans-Golgi, but not with the cis-Golgi apparatus (By similarity). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images