



S26A1 rabbit pAb

Catalog No	YP-Ab-09581
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC26A1 SAT1
Immunogen	Synthesized peptide derived from human S26A1 AA range: 502-552
Specificity	This antibody detects endogenous levels of human S26A1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S26A1;SLC26A1 SAT1;Sulfate anion transporter 1 (SAT-1) (Solute carrier family 26 member 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Basolateral cell membrane .
Tissue Specificity	Expressed most abundantly in the kidney and liver, with lower levels in the pancreas, testis, brain, small intestine, colon, and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
