



ACTZ rabbit pAb

Catalog No	YP-Ab-09585
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACTR1A CTRN1
Immunogen	Synthesized peptide derived from human ACTZ AA range: 260-310
Specificity	This antibody detects endogenous levels of human ACTZ
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACTZ;ACTR1A CTRN1;Alpha-centractin (Centractin) (ARP1) (Actin-RPV) (Centrosome-associated actin homolog)
Observed Band	41kD
Calculated Molecular Weight	41kD
Cell Pathway	Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cell cortex .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
