



KR103 rabbit pAb

Catalog No	YP-Ab-09587
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KRTAP10-3 KAP10.3 KAP18-3 KRTAP10.3 KRTAP18-3 KRTAP18.3
Immunogen	Synthesized peptide derived from human KR103 AA range: 27-77
Specificity	This antibody detects endogenous levels of human KR103
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KR103;KRTAP10-3 KAP10.3 KAP18-3 KRTAP10.3 KRTAP18-3 KRTAP18.3; Keratin-associated protein 10-3 (High sulfur keratin-associated protein 10.3) (Keratin-associated protein 10.3) (Keratin-associated protein 18-3) (Keratin-associated protein 18.3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Restricted to a narrow region of the hair fiber cuticle, lying approximately 20 cell layers above the apex of the dermal papilla of the hair root; not detected in any other tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images