



CH014 rabbit pAb

Catalog No	YP-Ab-09588
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LINC00208 C8orf14 NCRNA00208
Immunogen	Synthesized peptide derived from human CH014 AA range: 10-60
Specificity	This antibody detects endogenous levels of human CH014
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CH014;LINC00208 C8orf14 NCRNA00208;Putative uncharacterized protein encoded by LINC00208
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

