



MANBA rabbit pAb

Catalog No	YP-Ab-09590
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MANBA MANB1
Immunogen	Synthesized peptide derived from human MANBA AA range: 166-216
Specificity	This antibody detects endogenous levels of human MANBA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Beta-mannosidase;GC4M13771;Lysosomal beta A mannosidase;Manba; MANBA_HUMAN;Mannanase;Mannase;MANBA MANB1;Beta-mannosidase (EC 3.2.1.25) (Lysosomal beta A mannosidase) (Mannanase) (Mannase)
Observed Band	
Calculated Molecular Weight	11kD
Cell Pathway	Lysosome .
Tissue Specificity	Ubiquitous. Detected in pancreas, kidney and placenta, ands at lower levels in liver, lung, brain, heart and muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images